



Kuehne
LABORATORY SERVICES

Ion Chromatography (IC) as the Right Tool

Learn about where to start with determining IC as the right tool, how it can replace more tedious methods, and what you'll need to get started.

David Miyamoto

Lab Technical Manager

Agenda

- Industries: **Chemical Manufacturing**, Pharmaceuticals, Academia / Research
- What do you Value (CTQs, Tolerance Levels)
- HPLC-IC
- Method Development
- Problem Solving / Continuous Improvements



Agenda

- 1 _____
- 2 _____
- 3 _____
- 4 _____
- 5 _____
- 6 _____



Industries

- Chemical Manufacturing
- Academia / Research
- Pharmaceuticals



Industries

| | Chemical Manufacturing | Pharmaceutical | Academia / Research |
|-------------|------------------------|------------------------------|------------------------|
| Focus | Robustness | Regulatory Compliance | Flexibility |
| Throughput | Moderate-High | Moderate-High | Low-Moderate |
| Sample Type | Known | Known, with more regulations | Diverse, often unknown |




Determining Value

Determine

- What the parameters are
- What are the specifications for the parameters
- The tolerance levels of those specifications
- The cost to value. Both direct and indirect



Critical to Quality (CTQ) Example



Sulfate Analysis | Indirect Cost Example

| | Gravimetric Method | | IC | |
|-------------|--|------------------|---|---|
| Sample Prep | Carefully reduce pH to <4. Heat until a boil. Add BaCl ₂ . Continue to let it boil. | 45 - 60 Minutes | Turn Instrument on. Let it equilibrate | 30 - 45 Minutes |
| | Let it cool. | 30 Minutes | Make a 1000x dilution. Place sample in Vial.* | 1 - 2 Minutes |
| | Weigh out a 0.47μ filter paper and weighing dish. Prepare the vacuum filter. | 5 - 10 Minutes | Prepare the Sample Table and put sample on autosampler. | 1 - 2 Minutes |
| Analysis | Using the filter vacuum, collect the precipitate. | 5 - 10 Minutes | Run Calibration* Check | 30 Minutes |
| | Dry Sample | 60 - 120 Minutes | Run Sample* | 30 Minutes |
| | Weight out Sample, Dish, Filer Paper. | 1 - 2 Minutes | | |
| Total Time | | ~ 230 Minutes | | ~ 110 Minutes (35 Minutes subsequent Runs) |



Ion Chromatograph

Where to Start

Contact the manufacturer. Tell them what you are analyzing for, at what quantities and in what type of sample. Next you will want to tell them how many samples a day you will be analyzing.

A Chlor Alkali Example

- Manufacturer Metrohm
- Ions of Interest Bromide, Chlorate, Chloride, Sulfate, Perchlorate
- Sample Type(s) Brine (NaCl), Sodium Hydroxide, Sodium Hypochlorite (Bleach)
- Range 0.2 ppm – 300 ppm



Metrohm IC 940

Extension 942

Adds another pump. Able to add another mobile phase.

858 Professional Sample Processor

Metrosep A Supp 19 Column + Metrosep A Supp 19 Guard Column



Methods | Example

Three (3) Main Methods

- Routine Method
- High Salt – Routine Method
- Perchlorate

Prioritized

Simplicity, easy switching between methods

What to Consider


- Injection Size
- Temperature
- Flow Rate
- Mobile Phase
- Column (already discussed)




Methods | Value | SIMPLICITY

| | Routine Method | Routine Method - High Cl | Perchlorate |
|--------------|---------------------------|---------------------------|--|
| Column | Metrosep A Supp 19 Column | Metrosep A Supp 19 Column | Metrosep A Supp 19 Column |
| Mobile Phase | Isocratic 3.6 mM | Isocratic 3.6 mM | Gradient (3 Steps) 1. 3.6 mM 2. 14.4 mM 3. 3.6 mM |
| Temperature | 30°C | 30°C | 30°C |
| Suppression | Normal | Normal | 1 Additional "Step" |

Switching Between Methods – Routine Method



Switching Between Methods – Perchlorate





Clear Results

The data is only good if everyone can interpret it. Prioritize communicating the species and the units.

Our results need to be Sodium Chloride (NaCl) in grams per liter (g/L)



Results Example

Cl = 190 ppm | How to get results to NaCl in g/L

Change ppm to mg/L

190 mg/L

Multiply by any Dilution Factors

In this case we diluted the original sample 1 : 1000 with DI water

190 mg/L of Cl * 1000 = 190,000 mg/L of Cl

More Unit Conversion

190,000 mg/L = 190 mg/L of Cl

One More Factor

MM of Cl = 35.45 g/mol

MM of NaCl = 58.44 g/mol

58.44 / 35.45 = 1.65

Get the Result

190 * 1.65 = 313.5 g/L of NaCl



Continuous Improvements

What to Change?


Flowrate

Mobile Phase


Sample Preparation

- Concentrations
- Additives?

Temperature



Mobile Phase | Changing Concentrations



5.4.3 Na_2CO_3 variation at constant NaHCO_3

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: -

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS

Temperature: 30 °C

Loop: 20 μL


Flow rate: 0.7 mL/min

Eluent:

- A) 0.25 mmol/L NaHCO_3 , 4.0 mmol/L Na_2CO_3
- B) 0.25 mmol/L NaHCO_3 , 6.0 mmol/L Na_2CO_3
- C) 0.25 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3
- D) 0.25 mmol/L NaHCO_3 , 10.0 mmol/L Na_2CO_3
- E) 0.25 mmol/L NaHCO_3 , 12.0 mmol/L Na_2CO_3

| Metrosep A Supp 19 - 150/4.0 | | mg/L |
|------------------------------|-----------|------|
| 1 | Fluoride | 10 |
| 2 | Chloride | 10 |
| 3 | Nitrite | 10 |
| 4 | Bromide | 10 |
| 5 | Nitrate | 10 |
| 6 | phosphate | 10 |
| 7 | Sulfate | 10 |

Mobile Phase | Changing Concentrations




There is a strong relation between the sodium carbonate concentration and the retention times of the anions. But notice not all anions are affected the same.

This could be useful if time is an issue when analyzing for sulfate or nitrate.

| Metrosep A Supp 19 - 150/4.0 | | mg/L |
|------------------------------|-----------|------|
| 1 | Fluoride | 10 |
| 2 | Chloride | 10 |
| 3 | Nitrite | 10 |
| 4 | Bromide | 10 |
| 5 | Nitrate | 10 |
| 6 | phosphate | 10 |
| 7 | Sulfate | 10 |

Mobile Phase | Adding Acetone



5.5.1 Variation of the acetone concentration

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS

Temperature: 30 °C

Loop: 20 µL


Flow rate: 0.7 mL/min

Eluent:

- A) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 0% acetone
- B) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 10% acetone
- C) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 20% acetone

| Metrosep A Supp 19 - 150/4.0 | | mg/L |
|------------------------------|-------------------|------|
| 1 | Monochloroacetate | 10 |
| 2 | Monobromoacetate | 10 |
| 3 | Dichloroacetate | 10 |
| 4 | Dibromoacetate | 10 |
| 5 | Trichloroacetate | 10 |

Mobile Phase | Adding Acetone




Adding acetone (or other organic solvents) could be useful especially when dealing with organic acids, hydrocarbons, etc.

The eluent will also be made more stable against bacterial contamination.

| Metrosep A Supp 19 - 150/4.0 | | mg/L |
|------------------------------|-------------------|------|
| 1 | Monochloroacetate | 10 |
| 2 | Monobromoacetate | 10 |
| 3 | Dichloroacetate | 10 |
| 4 | Dibromoacetate | 10 |
| 5 | Trichloroacetate | 10 |

Change in Temperature



5.5.1 Variation of the acetone concentration

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS

Temperature: 30 °C

Loop: 20 μL


Flow rate: 0.7 mL/min

Eluent:
 A) 0.25 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3 , 0% acetone
 B) 0.25 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3 , 10% acetone
 C) 0.25 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3 , 20% acetone

Time [min]

| | Metrosep A Supp 19 - 150/4.0 | mg/L |
|---|------------------------------|------|
| 1 | Monochloroacetate | 10 |
| 2 | Monobromoacetate | 10 |
| 3 | Dichloroacetate | 10 |
| 4 | Dibromoacetate | 10 |
| 5 | Trichloroacetate | 10 |

Change in Temperature




In the monovalent anions (F, Cl, NO₂, Br, NO₃) retention times (RT) decrease with the increase of temperatures, while the multivalent anion's RT increases.

Also, the Multivalent peaks become sharper.

| Metrosep A Supp 19 - 150/4.0 | | mg/L |
|------------------------------|-------------------|------|
| 1 | Monochloroacetate | 10 |
| 2 | Monobromoacetate | 10 |
| 3 | Dichloroacetate | 10 |
| 4 | Dibromoacetate | 10 |
| 5 | Trichloroacetate | 10 |

Change in Flow Rate



5.3 Variation of the eluent flow rate

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS

Temperature: 30 °C


Loop: 20 μL

Flow rate: 0.3–0.9 mL/min

Eluent: 0.25 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3

| | Metrosep A Supp 19 - 150/4.0 | mg/L |
|---|------------------------------|------|
| 1 | Fluoride | 10 |
| 2 | Chloride | 10 |
| 3 | Nitrite | 10 |
| 4 | Bromide | 10 |
| 5 | Nitrate | 10 |
| 6 | phosphate | 10 |
| 7 | Sulfate | 10 |

Change in Flow Rate




Increasing the flowrate will decrease the retention times for all the anions but at the expense of the peak's area / height.

| Metrosep A Supp 19 - 150/4.0 | | mg/L |
|------------------------------|-----------|------|
| 1 | Fluoride | 10 |
| 2 | Chloride | 10 |
| 3 | Nitrite | 10 |
| 4 | Bromide | 10 |
| 5 | Nitrate | 10 |
| 6 | phosphate | 10 |
| 7 | Sulfate | 10 |

Sample Prep


Top Example

Chromatogram of Sodium Hydroxide sample without pretreatment



Bottom Example

Results after passing the Sodium Hydroxide sample through a cartridge





Kuehne
LABORATORY SERVICES

Contact Us

Phone: 973 589-0700

Email: labservices@kuehnecompany.com

Website: www.kuehnecompany.com/laboratory-services

Address: Corporate Headquarters

86 N Hackensack Avenue

Kearny NJ 07032

